

# Pillows

# SleepAngel®


## Safe and Sustainable Choice

Pillow supports the head at a comfortable level and ensures proper spine alignment. It is positioned closest to the sleeper's eyes, nose, mouth, and ears.

**Providing a safe and supportive pillow is, therefore, a must.**

SleepAngel pillows with PneumaPure® filter technology deliver this by providing your patients with premium comfort and unbeatable hygiene, while also reducing the environmental footprint through long-term usability.

- **Premium comfort** - SleepAngel pillows are used in hospitals and care homes around the world.
- **Infection prevention and control** - hermetically sealed positioners are ventilated through PneumaPure® filter and are a barrier to liquids, pathogens, and allergens.
- **Sustainable choice** - products can be used safely longer than traditional pillows and do not require laundry or dry cleaning.



PneumaPure® Filter technology



### No stains, no smells

always as new pillow

### No laundry or dry cleaning

wipe clean, disinfect, no need for pillow protectors

### PrimaLoft® Bio™ PRIMALOFT. BIO™

sustainable biodegradable synthetic fibre filling option

“

*Sustainable choice for dignified bed environment. SleepAngel combines premium comfort with unique infection prevention and control functionalities.*

# PneumaPure® filter technology

# SleepAngel®

PneumaPure® filter is a highly specialised nano-porous filter composite (i.e., composed of extremely small pores) that functions as an effective barrier to pathogens including bacteria, virus, and fungus as well as to the ingress of liquid, while remaining highly air-permeable.



## Barrier components

**PureShield** fabric protection: durable, soft touch and high performance textile, that blocks out contaminants. Vapour permeable and waterproof;

**SafeWeld** seams: heat sealed seams provide complete barrier – no stitching holes;

**PneumaPure®** filter: breathable mechanical filter allows clean air flow into the pillow, but blocks out pathogens and allergens.



Certified as Class I Medical Device in the EU  
Registered Class I Medical Device within the FDA  
Approved and registered with the SFDA

PneumaPure® Filter technology is patent protected globally, independently tested and proven by Airmid Healthgroup. Product range includes pillows, mattresses, positioners, blankets, top mattresses, and other soft surface products.

## Inspection and cleaning

SleepAngel Barrier Bedding alongside an audit and cleaning protocol will secure safe environment for the patient.

- **Audit** the product regularly – visually inspect the product's seams, fabric and filters for damage or holes. A product with holes or damaged fabric or filter must be disposed of.
- **Clean** the product at least once a week, after each patient discharge or when soiled. Clean with mild detergent and water solution.
- **Disinfect** the pillow if needed.

Detailed cleaning instructions accompany the product.



## Filling materials

We use innovative and proven materials and work with practitioners globally to find the best solutions.



## Contact

SleepAngel Medical is  
produced by Gabriel Scientific

 [www.SleepAngel-Medical.com](http://www.SleepAngel-Medical.com)

 [info@SleepAngel-Medical.com](mailto:info@SleepAngel-Medical.com)